

Light Sources

Product	Recommended Illumination Width	Bulb Type and Wavelength Range	Power	Part Number
<i>Dual 4-bulb (20W) line lights</i>	5-6"	Halogen 400-2500 nm	160 W - lights 300 W - power supply	MRC-920-029
<i>Dual 4-bulb (35W) line lights</i>	8"	Halogen 400-2500 nm	280 W - lights 300 W - power supply	MRC-920-031
<i>Dual 4-bulb (35W) adjustable line lights</i>	8"	Halogen 400-2500 nm	280 W - lights 960 W - power supply	MRC-920-031-II
<i>Dual 8-bulb (35W) line lights</i>	17-18"	Halogen 400-2500 nm	560 W - lights 960 W - power supply	MRC-920-027
<i>Single 16-bulb (35W) line light</i>	36"	Halogen 400-2500 nm	560 W - lights 960 W - power supply	MRC-920-030
<i>Dual 12-bulb (35W) line lights</i>	28"	Halogen 400-2500 nm	840 W - lights 1120 W - power supply	MRC-920-033
<i>Dual 13-bulb (35W) line lights</i>	31"	Halogen 400-2500 nm	910 W - lights 1120 W - power supply	MRC-920-034
<i>Dual LED array</i>	8"	LED 400-800 nm	84 W - lights	MRC-920-036
<i>Diffuse Transmission 4-bulb (20W) line light</i>	6"	Halogen 400-2500 nm	80 W - lights 300 W - power supply	MRC-920-029-DT
<i>Diffuse Transmission 4-bulb (35W) line light</i>	8"	Halogen 400-2500 nm	140 W - lights 300 W - power supply	MRC-920-031-DT
<i>Diffuse Transmission 8-bulb (35W) line light</i>	17-18"	Halogen 400-2500 nm	280 W - lights 960 W - power supply	MRC-920-027-DT
<i>Diffuse Illuminator</i>	28"	Halogen 400-2500 nm	315 W - lights 960 W - power supply	MRC-920-037-DT